# **Shipping Instructions for your Litholink Test**

## **Three Easy Steps**

Use the three steps below to schedule your one-time pickup:

#### Step 1

Call FedEx at 1.800.463.3339. You will be asked a few questions.

FedEx: "In a few words, please tell me what you're calling about..."

Your Answer: FedEx Express® return pickup

FedEx: "Is the word "stamp" written anywhere on your return label?"

Your Answer: No

FedEx: "Please say or enter the tracking number on your return label."

Your Answer: Slowly speak, or enter the tracking number on your return box

**FedEx:** "Before we move on to the pick-up address, the company name I have on

file is Litholink/Labcorp. Is the company name at the pick-up location?"

Your Answer: No

When you say "no" it will then prompt you to enter YOUR pick-up address.

#### Step 2

To be completed on the phone menu, provide YOUR:

- Phone number at pick-up location
- · Pickup location full address (include apartment/unit/suite)
- Number of packages to be picked up (1), and approximate weight (1lb)
- The day of the pick-up, when will it be ready, and how late we can pick up

### Step 3

Once your pick-up address is confirmed, it will then provide you with a pick-up confirmation number. Note this confirmation pickup number provided.

Also, make note of the 12-digit tracking number on your return box.



If the phone menu does transfer you to a live agent for whatever reason, please remind the agent that you would like to schedule a <u>FedEx Express® return pickup</u> and you already have the label and tracking number. No account numbers required.

